



With PRISMA 2, you simply hear better because ...

- The perfect combination of innovation and technology that has proven itself many times over.
- State-of-the-art PRISMA 2 electronics effectively suppresses disruptive noise.
- A highly intelligent hearing instrument that automatically adjusts to accommodate your life.
- PRISMA 2 means fewer restrictions on your ability to communicate – not even when you are using a cell or cordless phone.
- Improved energy management means PRISMA 2 uses even less power.
- Working together with you, your hearing care professional can individually custom-fit PRISMA 2 to your hearing loss.

Recommended by

Siemens
 Audiologische Technik GmbH
 Gebbertstrasse 125
 91058 Erlangen
 Germany
 Tel. +49 (9131) 308-0

The information in this document contains general descriptions of the technical options available, which do not always have to be present in individual cases and are subject to change without prior notice. The required features should therefore be specified in each individual case at the time of conclusion of the respective contract.

www.hearing-siemens.com

Order No. A91100-MS100-B947-01-7600
 Printed in Germany
 1043C6934 58 01035

PRISMA 2 will keep you right in the thick of things despite your hearing loss. We designed every detail of PRISMA 2 to maximise your comfort, which allows you to enjoy life to the fullest.

Small, Smart and Sturdy



You certainly do not want your hearing loss to be the first thing people notice. That is why we have made the PRISMA 2 Behind-The-Ear (BTE) hearing instruments so small and inconspicuous. And the tiny In-The-Ear (ITE) units are hardly noticeable. Just ask your hearing care professional whether you, too, can enjoy the benefits of a powerful yet inconspicuous In-The-Ear hearing instrument.

Here are a few examples from our range of products. PRISMA 2 is available in a number of different colors and versions.



Micro-CIC PRISMA 2 CS PRISMA 2 PRISMA 2 SP+

Always know what is going on

Being able to follow conversations better is not your only advantage with PRISMA 2. The system also keeps you informed about its status by emitting soft signals to tell you when you have switched programs, for example, or to let you know when it is time to change the batteries again. Plus, the new PRISMA 2 now features even lower power consumption.

With PRISMA 2, annoying whistling tones when you insert your hearing instrument are now a thing of the past. A special programmable soft-start function waits a few seconds after the unit is turned on before activating the amplifier, thus eliminating the unpleasant feedback noise.



The Right Move for Better Hearing: PRISMA 2



Siemens Hearing Solutions

SIEMENS

For years, PRISMA™ hearing instruments from Siemens have always been one move ahead. Today, they continue to be the global benchmark against which other directional microphone hearing instruments are measured.

Ahead of its
Time



PRISMA 2 takes a major step forward which means more hearing comfort for you than ever before! With an entire series of improvements, we have made sure that PRISMA 2 gives you everything you need for great hearing today and in the future.



Bonus Points for Better Hearing

With hundreds of thousands sold worldwide, the PRISMA message could not be clearer – Siemens' 125 years of experience in the field of audiology stands for solid expertise and has helped us to win the trust of many customers.

The best just keeps getting better

Now, PRISMA 2 has arrived on the scene. Our successful hearing instrument has improved even more, and offers you significantly greater listening, fitting and wearing comfort. With an entire line of Behind-The-Ear (BTE) and In-The-Ear (ITE) models, PRISMA 2 offers exactly the right solution for your hearing loss. In short, PRISMA 2 is now better than ever!

First impressions count, so we want you to feel good about PRISMA 2 right from the start. And we make sure you will, with our especially fast and comfortable fitting process that you actively participate in.

Better understanding with TwinMic

Our directional microphone technology has always been pioneering, but we have kept right on working to make it even better. The directional microphone systems (TwinMic™) in PRISMA 2 can switch back and forth between omnidirectional and directional mode to accommodate your hearing situation. This provides even better performance in noisy situations.



You can rely on our expertise to help you hear better

Our speech-sensitive signal processing improves comprehension, even in critical situations, and dampens disruptive background noise. And, best of all, you do not even need to adjust the volume, because PRISMA 2 automatically sets itself to the perfect level!

PRISMA 2 also provides you with excellent protection against annoying feedback. The optimised feedback management system excels at filtering out the whistling noises that occur in critical situations.